

2009 SCOOTER COMPARO!

MOBILITY

M A N A G E M E N T

February 2009 • Vol. 8 No. 2 Driving Today's Mobility & Rehab Markets

NIMEDA

New Products



Removable Seat Cover

Retractable fabric temporarily covers automotive seats to offer upholstery protection from dirt, grime and moisture. Fabric can be easily wiped clean; entire system can be removed quickly and stored conveniently under the seat or in the trunk. Fabric is designed to be water-resistant and is treated with an antibacterial agent to minimize odors.

Save-A-Seat
www.armarko.com



Superarm 2000E Lift

Compact lifting-arm system does not use a platform; wheelchairs are lifted out of the vehicle and deposited right beside it. Superarm does not block windows or doors, and can be reinstalled in subsequent vehicles. The lift can be operated by the user while seated in the wheelchair, or by an attendant. Superarm has a 600-lb. weight capacity and can fit into minivans and full-size vans.

Handicaps, Inc.
(800) 782-4335
www.handicapsinc.com

Save-A-Seat gains interest from people managing a wide range of disabilities.